

¹³C HeliFANplus and accessories



¹³C-breath test analyser
HeliFANplus

To complete the FANci device set the **HeliFANplus** has been developed. It is a compact and very reasonable ¹³C-infrared analyser with a high sample throughput. It is especially designed for performing the Helicobacter pylori breath test at the clinics and at consultants' office. Together with this device you obtain the capable and easy to operate **FANci**-software, including more than 20 additional breath tests. Due to the availability of all common interfaces (USB, CAN, COM, LPT) you can connect **HeliFANplus** to your PC or laptop without any problems.

Order

F301-HF-20




Accessories

Breath bags and mouthpieces for sampling, fitting the ¹³C-breath test analysers

Order

- | | |
|-----------------------------------|-------------|
| breath bag 0.3 l | F201-VP-05a |
| double breath bag | F201-VP-05b |
| breath bag 1.5 l | F201-VP-05c |
| mouthpiece (100 pieces) | F201-VP-06 |
| mouthpiece with valve (25 pieces) | F201-VP-11 |




Breathing masks

Breathing mask with check valve and bag adaptor for sampling with babies, disabled persons and for intensive care

Order

- | | |
|---------------------------|------------|
| breathing mask for babies | F201-VP-08 |
| breathing mask, medium | F201-VP-09 |
| breathing mask, large | F201-VP-10 |



Feel free to consult us about getting substances labelled with stable isotopes and to ask for further application of our ¹³C-analysers.