



Medtronic

BRAVO™ pH MONITORING SYSTEM

Catheter-Free pH Testing



*"The Bravo was invisible.
No one knew I was having
the test done... I drove and
we went shopping and had
dinner that night."*

Genevieve



Medtronic Gastroenterology
4000 Lexington Ave. N.
Shoreview, MN 55126-9866 USA
Telephone: (763) 514-9700
FAX: (763) 514-9745
Toll-free: (800) 227-3191
www.medtronic.com

UC 200300235
©Medtronic, Inc. 2002
All Rights Reserved

Medtronic Gastroenterology. One Source. Many Solutions.

Physicians and Patients Applaud BRAVO™ Catheter-Free pH Testing.

- Bravo can be placed during endoscopy
- Allows 24-hour or 48-hour pH monitoring to record more clinical data
- Patients maintain normal diet and routine activities
- Minimizes throat and nasal discomfort
- Transmits data to pager-sized receiver
- Uploads easily to pH analysis software

Easy to place during endoscopy, Bravo puts pH testing in the hands of every GI physician.



A Simple, Efficient Procedure.

Miniature Bravo pH capsule the size of a gelcap is temporarily attached to the wall of the esophagus.

Bravo capsule transmits pH data via radio frequency to a pager-sized receiver worn by the patient.

Data is uploaded to Polygram.Net Analysis Software, providing a complete report for patient diagnosis.

Bravo capsule naturally sloughs off the esophagus after a few days and passes through the digestive tract.



BRAVO lets patients go about their daily routines without interruption.



Revolutionizing pH Testing.

The catheter-free Bravo™ pH Monitoring System is revolutionizing the way pH testing is done. Bravo allows patients to maintain their regular diet and activities during pH testing. And because it is catheter-free, the added convenience for patients allows more clinical information to be recorded for physicians.

Easy to place during endoscopy, Bravo puts pH testing in the hands of every GI physician. With 24-hour or 48-hour monitoring and powerful diagnostic software, Bravo provides the advanced, yet simple-to-use tool physicians need for accurately diagnosing GERD.

One Source. Many Solutions.

Medtronic Gastroenterology is your one source for the newest and most innovative diagnostic and therapeutic solutions. As the world's leader in pH and motility, Medtronic offers the industry's broadest range of advanced products and treatment options.

From GI diseases including GERD and gastroparesis to pelvic floor disorders Medtronic offers the state-of-the-art technologies and comprehensive clinical training you need to improve the lives of your patients.

BRAVO™

Note: Reference product directions for use including indications, contraindications, warnings and precautions.
Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.



Medtronic
When Life Depends on Medical Technology

www.medtronic.com