



Medtronic

SLIMLINE™

Single-Use pH Catheter


SLIMMER
SMOOTHER
SOFTER

Slimline single-use pH catheter provides comfort with its slim body and frosted tubing.

Slimline offers:

Enhanced Patient Comfort and Safety:

- Slimmer: 1.8 mm tubing body to provide patient comfort
- Smoother: Frosted tubing creates slick surface to ease intubation
- Softer: Soft tubing body that softens further with body temperature for patient comfort



Moisture pockets on the catheter surface trap liquid molecules to increase wettability and facilitates intubation



Medtronic Gastroenterology. One Source. Many Solutions.



SLIMLINE™ Single-Use pH Catheter



Additional Benefits Include:

- Firm tip for more accurate positioning that minimizes curling during intubation
- Single-use for convenience that may reduce the risk of infection
- Clear, durable markings and numbers to facilitate accurate catheter placement
- Validated for use with Synectic and Medtronic Digitrapper (may be compatible with other pH recorders)
- Durable 8-pin modular connector for a more secure fit
- Time-proven accuracy with Medtronic antimony sensor
- Single channel and several dual channel configurations



Catalog #	Sensor(s)	20 cm	15 cm	10 cm	5 cm	0 cm	Sensor Spacing	
9012P3001	1	[Bar with 1 sensor at 0 cm]					1	None
9012P3021	2	[Bar with 2 sensors at 5 cm and 0 cm]				1	5 cm	
9012P3031	2	[Bar with 2 sensors at 10 cm and 0 cm]			1	10 cm		
9012P3041	2	[Bar with 2 sensors at 15 cm and 0 cm]		1	15 cm			
9012P3051	2	[Bar with 2 sensors at 18 cm and 0 cm]	1	18 cm				
9012P3061	2	1	21 cm					

Other models available including above configurations with added LES Indicator
(Sensor placement may vary by product)



Medtronic

When Life Depends on Medical Technology

Medtronic Gastroenterology

4000 Lexington Ave. N. MS X180
Shoreview, MN 55126-2983
USA

Toll-free: (866) 742-7286

Fax: (800) 930-9361

e-mail: slimline.ph@medtronic.com

www.medtronic.com

Medtronic A/S

Tonsbakken 16-18
DK-2740 Skovlunde
Denmark

Phone: 45 44 57 90 00

Fax: 45 44 57 90 10

www.medtronic.com

UC200303950 EN N15702

© Medtronic, Inc. 2003

All Rights Reserved

Printed in USA